

BRAND NAME	DESCRIPTION	MANUFACTURER	PURCHASE SOURCE
EYE PROTECTION			
Efficacy Criteria: Completely cove	rs eyes to block splatter, tolerates routine disinfe	ection, comfortable.	
Face Shield For Operator			
Loupes Face Shield	Shield that accommodates loupes & headlamp	Ultra Light Optics	ultralightoptics.com
<b>Eye Covering for Operator or</b>	Patient		
Googles Eye Shields	Wrap around disposable plastic lens & frame	KaVo Kerr Group	Most Dental Distributors
Helpful Accessory	,	l	
Kimtech Lens Cleaning Station	Kit= Cleaning solution & non-scratching Kimwipes	Kimberly Clark	Most Dental Distributors
FACE MASKS			
Efficacy Criteria: Adapts well to fo	acial contours (peripheral fit), high filtration of ve	ery small particles, does i	not contact nose or lips.
Soft Pleated Masks			
Critical Cover PFL with Magic Arch	Only Level 3 mask with potential for 360° excellent facial fit on most faces (Headband 608 series) or (Earloop 615 series)	Alpha ProTech	800-749-1363 Alphaprotech.com
Isolator Plus N95 Mask	N95 certified respirator has same design as Critical Cover, but upgraded filtration & double headband	Crosstex	Most Dental Distributors
Preformed Cup Masks		T	
Triple Layer Molded Mask #2042	Pre-formed cup & headband with good fit	3M	Benco Dental
#1860 N95 Molded Mask	N95 certified pre-formed cup mask	3M	Benco Dental
GLOVES			
Efficacy Criteria: Low defects, god	od fit, good tactile sensitivity, resistance to tearin	ng, adequate cuff length.	
Xlim Nitrile Gloves	Nitrile powder-free exam gloves	Cranberry	henryschein.com
Nitrident Stretch+	Nitrile powder-free exam gloves	Am-Touch	amtouch.com
AloePRO Synthetic Exam Gloves	Textured <u>latex</u> powder-free exam gloves	Dash Medical	dashmedical.com
Le Soothe Sapphire	<u>Latex</u> powder free exam gloves	SmartPractice	smartpractice.com
PROTECTIVE CLOTHING THA	AT IS AUTOCLAVABLE & DISPOSABLE		
Efficacy Criteria: Made of polypro	pylene to tolerate repeated steam sterilization		
Knee Length Disposable Lab Coat	SMS fabric coats come in 12 colors	Sara Glove	866-664-7272 saraglove.com
Lab Jacket, Disposable	SMS fabric coats come in 30 colors	Sara Glove	866-664-7272
DDE TREATMENT DINICE			saraglove.com
PRE-TREATMENT RINSE			
	eduction of relevant microbes <u>within</u> clinically rel 100 ppm <u>molecular</u> iodine ready-to-use rinse		561-509-0205 Ext. 5
ioRinse RTU	(use two 30-second rinses)	ioTech International	lotechinternational.com
Peridex	0.12 chlorhexidine gluconate	3M	Most Dental Distributors
DENTAL UNIT WATERLINE	VALIDATION		
Efficacy Criteria: For validation of	heterotrophic organism counts in water emitted	d from unit syringes & ha	ndpieces.
Independent Lab that Verifie	s Dental Unit Water by Culturing Water	<b>Samples Submitted</b>	by Office
Standard Check D-R2A	Office orders online number of sample tubes wanted & receives kit with directions. & packaging for sending samples to lab for testing.	Sterisil	719-662-7200 sterisil.com
"Quick Look" Dental Unit Wa	ter Culture Performed by Staff In-Office		
Quick Pass In-Office Waterline Tester	Culture medium on paddle, shows results in about 72 hours.	ProEdge Dental Products	henryschein.com

AIR PURIFICATION DURING	ACTIVE TOOTH EXCAVATION				
	e high efficacy air purification <u>during active rota</u>	ry instrument use & requi	res minimal maintenance.		
ECM Phantom	2x4 foot fan filter unit provides HEPA pure air; ceiling hung over treatment area.	Gordon Cleanroom	888-315-1561 gordon-inc.com		
SURFACE DISINFECTANTS			,		
Efficacy Criteria: Verified fast kill o	of resistant non-enveloped virus in presence of f	resh human whole blood.			
ZaPro Surface Disinfectant	95% wt ethyl alcohol, 0.13 wt benzalkonium chloride, 0.005% wt zinc pyrithione (1 min contact time)	Cao Group	877-877-9778 caogroup.com		
BioSURF	70.5% wt ethyl alcohol, .02% chlorhexidine gluconate (3 min contact time)	Micrylium	877-777-3303 www.purelifedental.com		
Lysol Brand III IC Disinfectant Spray	58% wt ethyl alcohol, 0.1% benzalkonium chloride (10 min contact time)	Reckitt Benckister	Most Dental Distributors		
HAND ANTISEPTICS					
Efficacy Criteria: Rapid antimici	obial activity & non-irritating to skin.				
Hand Wash 4% Chlorhexidine	Formulation				
Hibiclens	4% chlorhexidine, 4% isopropyl alcohol	Mölnlycke Health Care	Most Dental Distributors		
Hand Rub Ethyl Alcohol Gel		•			
ZaPro Hand Sanitizer	94% wt ethyl alcohol gel .05% zinc pyrithione in cellulose. Dries on hands quickly.	Cao Group	877-877-9778 caogroup.com		
Purell Advanced	70% volume ethyl alcohol gel	GoJo	Most Dental Distributors		
INSTRUMENT CLEANING &	DISINFECTION				
Efficacy Criteria: Instruments clea whole blood & other human protein	aned without hand touching & disinfected with ans.	chemicals verified to kill in	n presence of fresh human		
<b>Instrument Washers for Clea</b>	ning Instruments Before Disinfection o	r Sterilization			
Dental Washer Disinfector	Under counter model #PG8581	Miele	mieleusa.com		
Hydrim L110W G4	Under counter or free-standing model	SciCan	scican.com		
3+% Glutaraldehyde Disinfec	tion				
Maxicide Plus	3.4% glutaraldehyde soak	Henry Schein	henryschein.com		
pdCARE Plus	3.4% glutaraldehyde soak	Patterson	pattersondental.com		
Non-Glutaraldehyde Disinfection Alternative for Individuals Sensitive to Glutaraldehyde					
Cidex OPA	0.55% o-phthalaldehyde	Advanced Sterilization	Most Dental Distributors		
STERILIZERS	-				
Efficacy Criteria: Provide steriliza	tion of instruments to meet FDA specifications,	with minimal damage to i	instruments.		
Steam Heat Sterilizer, Rapid					
StatIM G4 (2000, 5000)	2000: steam; small chamber (107.8 in <sup>3</sup> ) 5000: steam; large chamber (315 in <sup>3</sup> )	SciCan	scican.com		
Steam Heat Sterilizer, Conver	ntional				
EZ11Plus	Steam; large chamber (1,882 in³)	Tuttnauer	tuttnauerusa.com		
M11 UltraClave	Steam; large chamber (1,710 in³)	Midmark	midmark.com		
STERILIZER MONITORS					
Efficacy Criteria: Biological monit	or that verifies heat sterilization.				
In-office Biological Monitor					
3M Attest Mini Auto-reader	Sterilization monitor set that gives <u>24-minute</u> results. (incubator & biological indicators)	3M	henryschein.com		