



# One Composite For All Situations: Excellent Universal Products

## Test Results

**Materials and methods:** Twenty-one total products tested, each given an overall score based on the following criteria (*in descending order of assigned importance*): ease of clinical placement (*handling and marginal adaptation*); polish retention (*following extended abrasion by standard brush and paste*); cost; polish at placement (*ease of polish*); flexural strength; depth of light cure; radiopacity; and translucency for esthetics (*compared to enamel*). All products listed are available in both compule and syringe dispensing versions.

Product Company	Approx. Cost/mL	Shades Available	Translucency for Esthetics (compared to enamel)	Ease of Placement (handling and marginal adaptation)	Radiopacity (% Aluminum equivalency; Dentin: 96; Enamel: 176)	Depth of Cure*	Flexural Strength (MPa) (higher is better)	Polish at Placement (GU) † (higher is better)	Polish Retention (GU) ‡ (higher is better)	Overall Score	Overall Rating
Harmonize Kerr	\$53	30	Excellent	Excellent	256	4.9mm	112	Excellent: 85	Excellent: 80	97%	Excellent
TPH Spectra ST Dentsply Sirona	\$59	6 (full VITA range)	Excellent	Excellent §	327	5.5mm	102	Excellent-Good: 60	Excellent: 88	96%	Excellent
Filtek Supreme Ultra 3M	\$65	36	Excellent	Excellent	261	5.0mm	125	Excellent: 87	Excellent: 88	95%	Excellent
Estelite Sigma Quick Tokuyama Dental	\$45	20	Excellent	Excellent	189	3.7mm	75	Excellent: 87	Excellent: 71	94%	Excellent
Beautiful II LS Shofu	\$34	14	Excellent-Good	Excellent-Good	230	6.0mm	105	Excellent: 78	Excellent: 63	93%	Excellent
Clearfil Majesty ES-2 Kuraray	\$43	18	Excellent	Excellent	151	4.7mm	85	Excellent: 73	Excellent-Good: 59	91%	Excellent
Venus Pearl Heraeus	\$65	27	Excellent	Excellent	276	4.4mm	155	Excellent-Good: 52	Excellent: 74	89%	Excellent
Mosaic Ultradent	\$53	20	Excellent	Good	263	5.3mm	150	Excellent: 74	Excellent: 70	87%	Excellent
G-aenial Sculpt GC America	\$57	17	Excellent-Good	Good	238	5.3mm	97	Excellent: 76	Excellent: 67	86%	Excellent-Good
Simile Pentron	\$35	20	Excellent	Good	214	5.5mm	117	Excellent-Good: 55	Good: 41	84%	Excellent-Good
Evanescence Clinician's Choice	\$58	31	Excellent	Good	223	5.0mm	98	Excellent: 80	Excellent-Good: 58	84%	Excellent-Good
Nuance DenMat	\$41	9	Excellent	Good-Fair	141	4.5mm	89	Excellent: 74	Excellent: 67	83%	Excellent-Good
Aura SDI	\$60	12	Excellent-Good	Good	229	5.5mm	91	Excellent: 72	Excellent-Good: 61	82%	Excellent-Good
Exquisite Restoration Apex Dental	\$36	8	Excellent	Fair	218	4.0mm	109	Excellent: 80	Good: 49	82%	Excellent-Good
Renamel Nano Plus Cosmedent	\$50	17	Excellent	Fair	276	5.4mm	102	Excellent: 67	Excellent: 63	81%	Excellent-Good
N'Durance Septodont	\$51	18	Excellent	Fair	382	5.3mm	105	Excellent: 77	Excellent-Good: 61	81%	Excellent-Good
Brilliant EverGlow Coltene	\$54	12	Excellent	Fair	198	5.8mm	102	Excellent: 83	Excellent: 69	78%	Excellent-Good
Esthet-X HD Dentsply Sirona	\$58	31	Excellent	Good	273	3.2mm	101	Excellent: 66	Good: 45	76%	Excellent-Good
Tetric EvoCeram Ivoclar Vivadent	\$53	22	Excellent	Fair	341	3.7mm	95	Excellent: 70	Good: 47	74%	Excellent-Good
Herculite Ultra Kerr	\$67	29	Excellent	Good-Fair	261	2.5mm	116	Excellent: 74	Excellent: 63	71%	Good
GrandioSO Voco	\$65	17	Excellent	Fair	255	3.3mm	109	Excellent-Good: 50	Good-Fair: 35	62%	Good

\* Using Valo Grand curing light by Ultradent (in standard mode) for manufacturer-recommended times, under idealized test conditions in extracted tooth

† Gloss reading using standard glossmeter, following initial polishing with series of progressively finer polishing grit levels, followed by polishing paste

‡ Gloss reading using standard glossmeter, following initial polishing and 40,000 cycles of abrasion by standard toothbrush and toothpaste

§ HV (high-viscosity) and LV (low-viscosity) versions available, depending on clinician preference.